



AbbaDox Interface HL7 Specification (Version 2.3)

Version 2.3

AbbaDox, Inc
20803 Biscayne Boulevard Suite 320
Aventura, FL 33180
Phone toll free 877-500-8614
Fax: 954-484-8491

Confidentiality Agreement

© Copyright 2020 – 2021 AbbaDox, Incorporated

All rights reserved. All information contained in this document is confidential and proprietary to AbbaDox Inc and intended only for the recipient to whom this copy is addressed. The recipient shall neither copy, distribute, disclose nor use the information contained in it for any purpose other than for which this document has been made available nor permit anyone else to do the same. The recipient shall be held fully liable for any wrongful disclosure or use of any of the information contained in this document by him, or his officers or employees or anyone else to whom he makes the information available.

www.abbadox.com



@AbbaDox #AbbaDox

ORM

Table of Contents

ORM - Message Definition	3
HL7 ORM Message Definition	4
Segment Table Definitions	4
Message Delimiters.....	5
MSH - Message Header Segment.....	6
PID - Patient Identification Segment.....	7
PV1 – Patient Visit Segment.....	9
ORC – Common Order.....	111
OBR – Order Details	12
Sample Messages.....	145

ORM - Message Definition

HL7 ORM Message Definition

An Order Message will consist of these segments:

MSH	Message Header	One per message
PID	Patient Identification	One per message
PV1	Patient Visit	One per message
ORC	Common Order	One per message
OBR	Observation Request	One per message

Segment Table Definitions

Heading	Contents	Values
SEQ	HL7 Field Sequence	Begins with '1' for each segment.
Element Name	Encoding Characters	Defined by HL7
DT	HL7 Data Type	Defined by HL7
R	Field Required	'R' = Required to process message. 'C' = Conditionally Required 'O' / Empty = Optional 'N' = Not supported
Comment	Field Usage Comments	

Message Delimiters

Delimiter	Suggested Value	Encoding Character Position	Usage
Segment Terminator	<cr> hex 0D	-	Terminates a segment record.
Field Separator		-	Separates two adjacent data fields within a segment. It also separates the segment ID from the first data field in each segment.
Component Separator	^	1	Separates adjacent components of data fields where allowed.
Subcomponent Separator	&	4	Separates adjacent sub-components of data fields where allowed.
Repetition Separator	~	2	Separates multiple occurrences of field where allowed.
Escape Character	\	3	Escape character for use with any field represented by a ST, TX or FT data type, or for use with the data (fourth) component of the ED data type

MSH - Message Header Segment

SEQ	Element Name	DT	R	Comment
1	Field Separator	ST	R	Fill in with ' '.
2	Encoding Characters	ST	R	Fill in with '^~\&'.
3	Sending Application	HD	R	IDSSEND
4	Sending Facility	HD	O	
5	Receiving Application	ST	O	
6	Receiving Application	HD	O	
7	Date/Time of Message	TS	R	yyyymmddhhmmss
8	Security	ST	O	
9	Message Type	CM_CMG	R	ORM^O01
9.1	Message Type	ID	R	ORM
9.2	Trigger Event	ID	R	O01
10	Message Control ID	ST	R	Identifier that uniquely identifies the message
11	Processing ID	PT	R	
12	Processing ID	ID	R	P
13	Version ID	ID	R	2.3

PID - Patient Identification Segment

SEQ	Element Name	DT	R	Comment
1	Set ID - Patient ID	SI	R	1
2	Patient ID (External ID)	CX	O	
3	Patient ID (Internal ID)	CX	R	MRN
4	Alternate Patient ID	CX	O	
5	Patient Name	XPN	R	Family Name^Given Name^Middle Name
5.1	Family Name	ST		
5.2	Given Name	ST		
5.3	Middle Name	ST		
6	Mother's Maiden Name	ST		
7	Date/Time of Birth	TS	R	The birth date of the patient (yyyymmdd).
8	Sex	IS	R	M/F/U
10	Race	IS	O	ID ^ Name
11	Patient Address	XAD	R	Address^Address2^City^State^Zip
11.1	Street Address	ST		
11.2	Other Designation	ST		
11.3	City	ST		
11.4	State	ST		
11.5	Zip	ST		
13	Phone Number - Home	XTN	O	
13.1	Telephone Number	TN		.

13.2	Telecommunication Use Code	ID		PRN / NET
13.3	Telecommunication Equipment Type	ID		PH/ INTERNET / CP / EMG
13.4	Email Address	ST		
14	Phone Number - Business	XTN	O	
14.1	Telephone Number	ST		
14.2	Telecommunication Use Code	ID		WPN
14.3	Telecommunication Equipment Type	ID		PH
15	Primary Language	CE	O	Id ^ Name
15.1	Identifier	ID		
15.2	Text	ST		
16	Marital Status	IS	O	
19	SSN	ST	O	
22	Ethnic Group	IS	O	Id ^ Name

PV1 – Patient Visit Segment

SEQ	Element Name	DT	R	Comment
1	Set ID	ID	R	1
2	Patient Class	EI	O	
3	Assigned Patient Location	PL	R	Location Code ^^Location Name
3.1	Point Of Care	ID		Location Code
3.4	Facility	HD		Location Name
4	Admission Type	ID	O	
5	Preadmit Number	CX		
6	Prior Patient Location	PL		
7	Attending Doctor	XCN	O	Attending Physician Id ^Last Name^First Name
8	Referring Doctor	XCN	R	Referring Physician Id ^Last Name^First Name^Middle Name
9	Consulting Doctor	XCN	O	Consulting Physician Id ^Last Name^First Name
10	Hospital Service	ID		
14	Admit Source	ID	O	
15	Ambulatory Status	IS		
16	VIP Indicator	ID		
17	Admitting Doctor	XCN	O	Admitting Physician Id ^Last Name^First Name
18	Patient Type	ID	O	
19	Visit Number	CX	O	
20	Financial Class	FC	O	

39	Servicing Facility	ID		
40	Bed Status	IS		
41	Account Status	ID		
42	Pending Location	PL		
43	Prior Temporary Location	PL		
44	Admit Date/Time	TS	O	yyyymmddhhmmss
45	Discharge Date/Time	TS	O	yyyymmddhhmmss

ORC – Common Order

SEQ	Element Name	DT	R	Comment
1	Order Control	ID	R	NW = New Order XO = Change Order CA = Cancel Order
2	Placer Order Number	EI	R	Accession Number
3	Filler Order Number	EI	O	
4	Placer Group Number	EI		
5	Order Status	ST	R	IP/CM/CA/RE
6	Response Flag	ST		
7	Quantity Timing	TQ		
8	Parent	CM_PA		
9	Date Time Of Transaction	TS		
10	Entered By	XCN		
11	Verified By	XCN		
12	Ordering Provider	XCN		
13	Enterers Location	PL		
13.1	Point Of Care	ID		

OBR – Order Details

SEQ	Element Name	DT	R	Comment
1	Set ID	SI	R	1
2	Placer Order Number	EI	R	Accession Number
3	Filler Order Number	EI		
4	Universal Service ID	CE	R	Exam Code^Exam Description
4.1	Identifier	ID		
4.2	Text	ST		
5	Priority	ID	O	R = Routine, S = Stat, A = ASAP
6	Requested Date Time	TS	O	yyyymmddhhmmss
7	Observation Date Time	TS	R	yyyymmddhhmmss Study Date and Time
8	Observation End Date Time	TS	O	yyyymmddhhmmss
9	Collection Volume	CQ		
10	Collector Identifier	XCN		
11	Specimen Action Code	ID		
12	Danger Code	CE		
13	Relevant Clinical Info	ST	O	
14	Specimen Received Date Time	TS		
15	Specimen Source	CM_SS		
16	Ordering Provider	XCN	O	Ordering Physician Mnemonic ^Last Name^First Name
16.1	ID Number	ST		
16.2	Family Name	ST		

16.3	Given Name	ST		
17	Order Callback Phone Number	XTN	O	
18	PlacerField1	ST	O	
19	PlacerField2	ST	O	
20	FillerField1	ST	O	
21	FillerField2	ST		
22	Results Rpt/Status Chng – Date/Time	TS		
23	ChargeToPractice	CM_CH		
24	Diagnostic Service Section ID	ID	R	Modality Code
25	ResultStatus	ST		
26	ParentResult	CM_PR		
27	Quantity Timing	TQ	R	^^^Study Date_Time^^Priority
27.1	Quantity	CQ		
27.2	Interval	CM_RI		
27.3	Duration	ST		
27.4	Start Date/Time	TS	R	yyyymmddhhmmss
27.5	End Date/Time	TS		
27.6	Priority	ST	R	R = Routine, S = STAT, A = ASAP
28	ResultCopiesTo	XCN	o	
29	ParentNumber	CM_PA		
30	TransportationMode	String		
31	Reason For Study	CE	o	^ Indication/Reason
31.1	Reason ID	ID	N	
31.2	Reason Text	ST	O	

Sample Message

```
MSH|^~\&|IDSSEND|||20140929132123||ORM^O01|2014092913212330886||2.3
PID|1||887654260||ZZTESTADULT^RegressionPGH^||19870101|F||ASPA^^ASPA^CA^81695||(888)888-
8888|||887654524|777777777|
PV1|1||1^^^Palmetto|||08888^GOLDSTEIN^MITCHELL|||08888^GOLDSTEIN^MITCHELL|||
|||||20140922083000|20140724235900|||||
ORC|NW|00349CT140000002||IP|||||
OBR|1|00349CT140000002||218057263^CT Bone Marrow Aspiration W/
Biopsy||20140924112500|||||09694^LOSADA
MD^ANGELES|||||CT||^20140924112500^^R||^Unexplained Anemia|||||
572^RAPPAPORT^KENNETH||0.2200|3|2||CT||^20140917171500^^R||^DX|||||
```