



AbbaDox Interface HL7 Specification (Version 2.3)

Version 2.3

AbbaDox, Inc
20803 Biscayne Boulevard Suite 320
Aventura, FL 33180
Phone toll free 877-500-8614
Fax: 954-484-8491

Confidentiality Agreement

© Copyright 2020 – 2021 AbbaDox, Incorporated

All rights reserved. All information contained in this document is confidential and proprietary to AbbaDox Inc and intended only for the recipient to whom this copy is addressed. The recipient shall neither copy, distribute, disclose nor use the information contained in it for any purpose other than for which this document has been made available nor permit anyone else to do the same. The recipient shall be held fully liable for any wrongful disclosure or use of any of the information contained in this document by him, or his officers or employees or anyone else to whom he makes the information available.

www.abbadox.com



@AbbaDox #AbbaDox

ORU

Table of Contents

ORU - Message Definition	3
MSH - Message Header Segment.....	4
PID - Patient Identification Segment.....	5
ORC – Common Order.....	7
OBR – Order Details	8
OBX –Observation	11
Sample Message.....	12

ORU - Message Definition

MSH - Message Header Segment

SEQ	Element Name	DT	R	Comment
1	Field Separator	ST	R	Fill in with ' '.
2	Encoding Characters	ST	R	Fill in with '^~\&'.
3	Sending Application	HD	R	IDSSEND
4	Sending Facility	HD	O	
5	Receiving Application	ST	O	
6	Receiving Application	HD	O	
7	Date/Time of Message	TS	R	yyyymmddhhmmss
8	Security	ST	O	
9	Message Type	CM_CMG	R	ORU^R01
9.1	Message Type	ID	R	
9.2	Trigger Event	ID	R	
10	Message Control ID	ST	R	Identifier that uniquely identifies the message
11	Processing ID	PT	R	
12	Processing ID	ID	R	P
13	Version ID	ID	R	2.3

PID - Patient Identification Segment

SEQ	Element Name	DT	R	Comments
1	Set ID	SI	R	1
2	Patient ID	CX	0	
3	Patient ID	CX	R	Patient Id
4	Alternate Patient ID	CX	0	
5	Patient Name	XP	R	Family Name^Given Name^Middle Name
5.1	Family Name	ST		
5.2	Given Name	ST		
5.3	Middle Name	ST		
6	Mother's Maiden Name	ST		
7	Date/Time of Birth	TS	R	The birth date of the patient (yyyymmdd).
8	Sex	IS	R	M/F/U
10	Race	IS	O	ID ^ Name
11	Patient Address	XAD	R	Address^Address2^City^State^Zip
11.1	Street Address	ST		
11.2	Other Designation	ST		
11.3	City	ST		
11.4	State	ST		
11.5	Zip	ST		
13	Phone Number - Home	XTN	O	
13.1	Telephone Number	TN		.
13.2	Telecommunication Use Code	ID		PRN / NET

13.3	Telecommunication Equipment Type	ID		PH/ INTERNET / CP / EMG
13.4	Email Address	ST		
14	Phone Number - Business	XTN	O	
14.1	Telephone Number	ST		
14.2	Telecommunication Use Code	ID		WPN
14.3	Telecommunication Equipment Type	ID		PH
15	Primary Language	CE	O	Id ^ Name
15.1	Identifier	ID		
15.2	Text	ST		
16	Marital Status	IS	O	
19	SSN	ST	O	
22	Ethnic Group	IS	O	Id ^ Name

ORC – Common Order

SEQ	Element Name	DT	R	Comment
1	Order Control	ID	R	XO
2	Placer Order Number	EI	R	Accession Number
3	Filler Order Number	EI	O	
4	Placer Group Number	EI		
5	Order Status	ST	R	RE
6	Response Flag	ST		
7	Quantity Timing	TQ		
8	Parent	CM_PA		
9	Date Time of Transaction	TS		
10	Entered By	XCN		
11	Verified By	XCN		
12	Ordering Provider	XCN		
13	Enterers Location	PL		
13.1	Point Of Care	ID		

OBR – Order Details

SEQ	Element Name	DT	R	Comment
1	Set ID	SI	R	1
2	Placer Order Number	EI	R	Accession Number
3	Filler Order Number	EI		
4	Universal Service ID	CE	R	Exam Code^Exam Description
4.1	Identifier	ID		
4.2	Text	ST		
5	Priority	ID	O	R = Routine, S = Stat, A = ASAP
6	Requested Date Time	TS	O	yyyymmddhhmmss
7	Observation Date Time	TS	O	yyyymmddhhmmss
8	Observation End Date Time	TS	O	yyyymmddhhmmss
9	Collection Volume	CQ		
10	Collector Identifier	XCN		
11	Specimen Action Code	ID		
12	Danger Code	CE		
13	Relevant Clinical Info	ST	O	
14	Specimen Received Date Time	TS		
15	Specimen Source	CM_SS		
16	Ordering Provider	XCN	O	Ordering Physician Mnemonic ^Last Name^First Name
16.1	ID Number	ST		
16.2	Family Name	ST		

16.3	Given Name	ST		
17	Order Callback Phone Number	XTN	O	
18	Placer Field 1	ST	O	
19	Placer Field 2	ST	O	
20	Filler Field 1	ST	O	
21	Filler Field 2	ST		
22	Results Rpt/Status Chng – Date/Time	TS	R	yyyymmddhhmmss
23	Charge To Practice	CM_CH		
24	Diagnostic Service Section ID	ID	R	Modality Code
25	Result Status	ST	R	F = Final, A = Addendum
26	Parent Result	CM_PR		
27	Quantity Timing	TQ	O	^^^Study Date/Time^^Priority
27.1	Quantity	CQ		
27.2	Interval	CM_RI		
27.3	Duration	ST		
27.4	Start Date/Time	TS	R	yyyymmddhhmmss
27.5	End Date/Time	TS		
27.6	Priority	ST	R	R = Routine, S = STAT, A = ASAP
28	Result Copies To	XCN	O	
29	Parent Number	CM_PA		
30	Transportation Mode	String		
31	Reason For Study	CE	O	^ Indication/Reason
31.1	Reason ID	ID	N	

31.2	Reason Text	ST	O	
32	Principal Results Interpreter	NDL	R	Rendering ID ^ Last Name ^ First Name
32.1	OPP Name	CN	R	
32.2	Start Date/Tie	TS	R	

OBX – Order Observation/Result

SEQ	Element Name	DT	R	Comments
1	Set ID	ID	R	
2	Value Type	ID	R	TX
3	Observation Identifier	CE	O	
4	Observation Sub-ID	ST	O	
5	Observation Value	FT	R	Each new line is in a new OBX segment Addendum to include Final result
6	Units	CE		
7	References Range	ST		
8	Abnormal Flags	ID		
9	Probability	NM		
10	Nature Of Abnormal Test	ID		
11	Observ Result Status	ID	O	F = Final A = Addendum
12	Date Last Obs Normal Values	TS		
13	User Defined Access Checks	ST		
14	Date/Time Of The Observation	TS	O	yyyymmddhhmm
15	Producers ID	CE		
16	Responsible Observer	XCN	O	Rendering Provider
16.1	ID Number	ST		
16.2	Family Name	ST		
16.2	Given Name	ST		

Sample Message

MSH|^~\&|IDSSSEND|||201107071443||ORU^R01|20110707144338239685|P|2.3
PID|1|887654260||ZZTESTADULT^ONE^||19870101|F|||ASPA^ASPA^CA^81695|||(888)888-8888|||887654524|777777777|
ORC|XO|00349CT140000002|||RE|||
OBR|1|00349CT140000002||218057263^CT abdomen pelvis with IV contrast|||20140924112500|||09694^LOSADA
MD^ANGELES|||201104042321||CT|F|^20140924112500^R|||^Unexplained Anemia|123456^Rad^Adam|||
OBX|1|TX|IOC0014835||Referring Physician: Dr. Smith|||F|||201104042321||123456^Rad^Adam|
OBX|2|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|3|TX|IOC0014835||Patient Name: ONE ZZTESTADULT|||F|||201104042321||123456^Rad^Adam|
OBX|4|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|5|TX|IOC0014835||CT abdomen pelvis with IV contrast April 4, 2011|||F|||201104042321||123456^Rad^Adam|
OBX|6|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|7|TX|IOC0014835||Images: 425|||F|||201104042321||123456^Rad^Adam|
OBX|8|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|9|TX|IOC0014835||Comparison: None|||F|||201104042321||123456^Rad^Adam|
OBX|10|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|11|TX|IOC0014835||History: Right lower quadrant abdominal pain, history of
stones|||F|||201104042321||123456^Rad^Adam|
OBX|12|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|13|TX|IOC0014835||Findings: There is no bowel obstruction, free air or free fluid in the abdomen. No localized inflammatory
changes are identified. The appendix is not identified. There are no findings to suggest acute appendicitis. The liver, spleen, pancreas,
adrenal glands, and kidneys appear unremarkable. Gallbladder is surgically absent. Abdominal aorta is normal in caliber. Urinary
bladder is mildly distended. There is dependent atelectasis in both lung bases.|||F|||201104042321||123456^Rad^Adam|
OBX|14|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|15|TX|IOC0014835||Impression: No evidence of acute intra-abdominal
process.|||F|||201104042321||123456^Rad^Adam|
OBX|16|TX|IOC0014835|||F|||201104042321||123456^Rad^Adam|
OBX|17|TX|IOC0014835||THIS DOCUMENT HAS BEEN ELECTRONICALLY SIGNED|||F|||201104042321||123456^Rad^Adam|